

Applied method (e.g. AOAC, DIN, EN, ISO, EPA, ASTM, §64, company sop, etc.)

DIN EN ISO method 14891, IDF 185, Milk and Milk Products – Determination of nitrogen content – routine method using combustion according to Dumas principle.

### Instruments

- |   |   |
|---|---|
| 1 | Analytical Balance (readability 0,1 mg or better)               |
| 2 | Homogenizer, e.g. Grindomix GM200 Knife Mill or a regular Knife |
| 3 | DUMATHERM N Pro, standard configuration                         |

### Gases and Consumables

- |   |   |
|---|---|
| 1 | Helium and Oxygen, bottle gas, min. quality grade 5.0             |
| 2 | Nitrogen or compressed air as bottle gas, min. quality grade 2.6. |
| 3 | DumaReact, Combustion Reactor, packed with catalyst, 14-0245      |
| 4 | DumaTube, Quartz glass for reactor, 14-0203                       |
| 5 | DumaFoil, Tin Foil for packing the samples, 14-0017               |
| 6 | DumaEDTA, Standard for Calibration, purity > 99 %, 14-0032        |

### Method Settings

Sample Weight	200 mg
Packing of the sample	Tin foil
Combustion Method	A 1,8 (400 ml O <sub>2</sub> / min, 1.8 ml O <sub>2</sub> / mg sample)
Protein Factor	6,38
Combustion Temperature	980 °C (alternatively 1030 °C)
Reduction temperature	650 °C (alternatively 750 °C)
Recommended Calibration Range	1 – 15 mg N absolute (measured with 10-150 mg EDTA)

### Homogenization / Preparation

The sample is taken from the fridge and directly homogenized with a regular kitchen grinder (type Moulinette). The mashed material is thoroughly mixed again with a spatula and taken with a spatula for weighing into the tin foil.

The material should be at room temperature ( $\leq 20$  °C) at most when it is weighed in.

Alternatively the sample is taken from the fridge and is knead at room temperature in the foil packaging to get a homogeneous sample material. The homogenized sample is then taken with a spatula for weighing.

Higher temperatures for combustion (1030 °C) and reduction (750 °C) will improve slightly the standard deviation. But as norm requirements are fulfilled already with the reduced temperatures, there is no need for the higher ones.

