**OPERATING MANUAL  
MODE D'EMPLOI  
INSTRUCCIONES DE OPERACIÓN**

T 814241 04/2014

TruLine ORP 50

**TruLine**

---

# TruLine

<b>Operating Manual .....</b>	<b>3</b>
Technical data.....	5
Commissioning, measuring, testing .....	6
Storage .....	8
Aging .....	8
Maintenance and cleaning .....	9
Activating platinum electrodes .....	10
Wear parts and accessories .....	11
<b>Mode d'emploi.....</b>	<b>13</b>
Caractéristiques techniques .....	15
Mise en service, mesure, vérification .....	16
Stockage .....	18
Vieillissement.....	18
Maintenance et nettoyage .....	19
Activation d'électrodes de platine.....	20
Pièces d'usure et accessoires.....	21
<b>Instrucciones de operación .....</b>	<b>23</b>
Especificaciones técnicas .....	25
Puesta en funcionamiento, medición, verificación .....	26
Almacenamiento .....	28
Envejecimiento .....	28
Mantenimiento y limpieza.....	29
Activar electrodos de platino .....	30
Piezas de desgaste y accesorios .....	31
<b>Contact Information.....</b>	<b>35</b>

**OPERATING MANUAL**

T 814241 04/2014

TruLine ORP 50

**TruLine**

ORP ELECTRODE

For the most recent version of the manual, please visit [www.ysi.com](http://www.ysi.com).

**Contact**

YSI  
1725 Brannum Lane  
Yellow Springs, OH 45387 USA  
Tel: +1 937-767-7241  
800-765-4974  
Email: [environmental@ysi.com](mailto:environmental@ysi.com)  
Internet: [www.ysi.com](http://www.ysi.com)

**Copyright**

© 2014 Xylem Inc.

## Technical data

### ***General data***

Model	Reference electrolyte	Electrode material and shape	Junction
TruLine ORP 50	3 mol/l KCl, Ag <sup>+</sup> -free	Platinum / Round piece	Ceramic

### ***Measurement and application characteristics***

Model	Allowed temperature range	Typical application
TruLine ORP 50	0 ... 100 °C	Laboratory

### ***Shaft dimensions, shaft material, electrical connection***

Model	Shaft			Electrical connection		
	Length [mm]	Ø [mm]	Material	Combination electrode connection	Meter connection	Cable length
TruLine ORP 50	120	12	Glass	Fixed cable	BNC	1 m

## Commissioning, measuring, testing

### Note


For ORP electrodes with platinum electrode, please observe the instructions in section ACTIVATING PLATINUM ELECTRODES.

### Commissioning

Prepare the combination electrode for measuring as follows:

- TruLine ORP 50: Open the refilling opening for the reference electrolyte solution.

**The refilling opening must always be open during measurement!**

- Remove the watering cap from the electrode tip. Possible salt deposits in the area of the watering cap do not affect the measuring characteristics and can easily be removed with deionized water.

### Note


Please keep the watering cap. It is required for storing the combination electrode. Always keep the watering cap clean.

- Connect the combination electrode to the meter.
- Measure with the combination electrode according to the operating manual of the meter and observe the following rules while doing so:

### General rules for measuring

- TruLine ORP 50: Make sure the refilling opening for the reference electrolyte solution is open.
- Avoid displacement of the sample solution from one measurement to the next by taking the following measures:
  - Briefly rinse the calibration beakers with the solution the beakers are to be filled with next.
  - Between measurements, rinse the combination electrode with the solution specified below. Alternatively, you can also rinse the combination electrode with deionized water and then carefully dab it dry.
- Immerse the combination electrode in the solution in a vertical or slightly tilted position.

- Make sure the immersion depth is correct. The junction must be completely submersed in the solution. The junction is in the area of the bottom end of the shaft (see arrow).

The level of the reference electrolyte must be at least 2 cm above the level of the solution.

#### Conversion to the standard hydrogen electrode potential

$$U_H = U_{\text{Meas}} + U_{\text{Ref}}$$

where:  $U_H$  = ORP voltage, relative to the standard hydrogen electrode

$U_{\text{Meas}}$  = measured ORP voltage

$U_{\text{Ref}}$  = voltage of the reference system relative to the standard hydrogen electrode

$U_{\text{Ref}}$  is temperature dependent and can be taken from the following table (also refer to DIN 38404-6):

T (°C)	U <sub>Ref</sub> [mV] TruLine ORP 50	T (°C)	U <sub>Ref</sub> [mV] TruLine ORP 50
0	+224	35	+200
5	+221	40	+196
10	+217	45	+192
15	+214	50	+188
20	+211	55	+184
25	+207	60	+180
30	+203		

## Storage

### During short measuring breaks

Immerse the electrode in the reference electrolyte with the refilling opening open.

Combination electrode	Reference electrolyte	YSI Model (see page 11)
TruLine ORP 50	3 mol/l KCl, Ag <sup>+</sup> -free	KCl-250 (250 ml)

Prior to the next measurement, briefly rinse the combination electrode with the test sample or deionized water.

### Overnight or longer

Insert the clean electrode into the watering cap filled with reference electrolyte and shut the refilling opening.

#### Note

During longer storage periods, salt deposits may develop on the watering cap. They do not affect the measuring characteristics and can easily be removed with deionized water when the combination electrode is put into operation again.

## Aging

Every ORP combination electrode undergoes a natural aging process. Extreme operating conditions can considerably shorten the lifetime of the combination electrode. These are:

- Strong acids or lyes, hydrofluoric acid, organic solvents, oils, fats, bromides, sulfides, iodides, proteins
- High temperatures
- Great changes in pH and temperature.

The warranty does not cover failure caused by measuring conditions and mechanical damage.


## Maintenance and cleaning

During operation, a small amount of reference electrolyte leaks through the junction from the combination electrode into the test sample. If the level of reference electrolyte becomes too low with time, refill it through the refilling opening.

### Refilling reference electrolyte

Refilling is very easy using a dropper bottle. To do so, proceed as follows:

- Cut off the tip of the dropping bottle at a right angle until the opening in the tip can be seen
- Open the refilling opening of the combination electrode
- Press the tip of the dropper bottle into the refilling opening while turning it slightly
- Pump several small quantities of the reference electrolyte into the stem using the dropper bottle
- Pull the dropper bottle out of the refilling opening while turning it slightly if necessary.

**Cleaning II** Remove water-soluble contamination by rinsing with deionized water. Remove other contamination as follows:

Contamination	Cleaning procedure
Fat and oil	Rinse with water containing household washing-up liquid
Lime and hydroxide deposits	Rinse with citric acid (10 % by weight)

**After cleaning** Rinse the combination electrode with deionized water.

## Activating platinum electrodes


### First-time activation during installation and as required

For the first-time activation use the activation powder from the SORT/RH reagent set (component of the SORT/RH reagent set). Immerse the moist (but not dripping) platinum electrode into the activation powder and rotate the electrode several times in the powder. Then remove the activation powder under flowing water with a soft brush (e.g. toothbrush).

### Activation during very long set-up times

When changing from oxidizing to reducing test solutions and vice versa this can result in set-up times that can take significantly more than an hour. In this case pretreatment (activation) of the platinum surface can shorten the set-up time. The type of pretreatment (reducing or oxidizing) is based on the pH value and the ORP voltage ( $U_H$ ) of the test solution where the latter must be estimated for the first measurement.

The type of pretreatment can then be determined using the following diagram where  $U_H$  is based on the normal hydrogen electrode:

<b>Oxidizing pretreatment</b>	Immerse the platinum electrode for two to three days in a sulfuric acid chlorine solution. Chlorine powder for producing the solution is included in the SORT/RH reagent set.  <u>Note:</u> The diaphragm must not be immersed in the chlorine solution!
<b>Reducing pretreatment</b>	When the electrode is ready for the test immerse it in the RH 28 ORP buffer solution and wait for a stable measured value.

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## Contact Information

### ***Ordering & Technical Support***

**Telephone:** (800) 897-4151  
(937) 767-7241  
Monday through Friday, 8:00 AM to 5:00 PM ET

**Fax:** (937) 767-1058

**Email:** [environmental@ysi.com](mailto:environmental@ysi.com)

**Mail:**  
YSI Incorporated  
1725 Brannum Lane  
Yellow Springs, OH 45387  
USA

**Internet:** [www.ysi.com](http://www.ysi.com)

When placing an order please have the following information available:

YSI account number (if available)	Name and Phone Number
Model number or brief description	Billing and shipping address
Quantity	Purchase Order or Credit Card

### ***Service Information***

YSI has authorized service centers throughout the United States and Internationally. For the nearest service center information, please visit [www.ysi.com](http://www.ysi.com) and click 'Support' or contact YSI Technical Support directly at 800-897-4151.

When returning a product for service, include the Product Return form with cleaning certification. The form must be completely filled out for an YSI Service Center to accept the instrument for service. The Product Return form may be downloaded at [www.ysi.com](http://www.ysi.com) and clicking on the 'Support' tab.

---

## What can Xylem do for you?

We're 12,500 people unified in a common purpose: creating innovative solutions to meet our world's water needs. Developing new technologies that will improve the way water is used, conserved, and re-used in the future is central to our work. We move, treat, analyze, and return water to the environment, and we help people use water efficiently, in their homes, buildings, factories and farms. In more than 150 countries, we have strong, long-standing relationships with customers who know us for our powerful combination of leading product brands and applications expertise, backed by a legacy of innovation.

For more information on how Xylem can help you, go to [www.xyleminc.com](http://www.xyleminc.com)

a **xylem** brand

YSI  
1725 Brannum Lane  
Yellow Springs, OH 45387  
Tel: +1 937-767-7241; 800-765-4974  
Fax: +1 937-767-1058  
Email: [environmental@ysi.com](mailto:environmental@ysi.com)  
Web: [www.yси.com](http://www.yси.com)

©Xylem Inc

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**Copyright** © 2014 Xylem Inc.