

# Automated Kappa Number Determination of Pulp

PULP & PAPER SERIES



## Introduction

The Kappa number describes the Lignin content of pulp which gives information about the bleaching process of the pulp. The longer the pulp was cooked and bleached the lower the Kappa number is. For the determination of the Kappa number a sample weight adjusted to the expected result is chopped and inserted into a beaker, then 400 ml of water is added. The sample/water mixture is stirred for 10 minutes. Afterwards 50 ml of potassium permanganate ( $\text{KMnO}_4$ , 0.02 mol/L) and 50 ml sulfuric acid ( $\text{H}_2\text{SO}_4$ , 2 mol/L) are added simultaneously. The mixture is stirred for another 10 minutes. The reaction is stopped by adding 10 ml potassium iodide (KI, 1 mol/l). The excess of  $\text{KMnO}_4$  reacts with Iodide to Iodine, which is titrated with sodium thiosulfate ( $\text{Na}_2\text{S}_2\text{O}_3$ , 0.2 mol/L).




a xylem brand

## Chemical Equations:


The weight is adjusted to a consumption of about 50% of the  $\text{MnO}_4^-$ .



The reaction is stopped by adding KI.



The formed Iodine is titrated with  $\text{Na}_2\text{S}_2\text{O}_3$ .



## Calculation of the Kappa number:

1. Calculation of the consumed volume  $V_a$  of  $\text{KMnO}_4$ :

$$V_a = \frac{V_1 - V_2)c}{0.1}$$

with

$V_1$  = Consumption of  $\text{Na}_2\text{S}_2\text{O}_3$  during blank titration

$V_2$  = Consumption of  $\text{Na}_2\text{S}_2\text{O}_3$  during sample titration

$c$  = Concentration of  $\text{Na}_2\text{S}_2\text{O}_3$

0.1 = numerical Factor, calculated from the molarity of the

permanganate and the stoichiometric factor  $f$  of the reaction:

$$C_{\text{KMnO}_4} \cdot f = 0.02 \cdot 5$$

2. Calculation of a correction factor  $d$ , which corrects the consumption of permanganate depending on  $V_a$  to a consumption of 50%.

$$d = 10^{0.00093(2V_a - 50)} \\ = e^{\ln(10) \cdot 0.00093(2V_a - 50)}$$

3. Calculation of the Kappa number corrected to 25 °C.

$$K = \frac{V_a \cdot d}{m} \cdot (1 + 0.013(25 - T))$$

with

$T$  = actual temperature, measured during the titration

$m$  = sample weight in g

Instrument	
TL 7000 with 20 ml interchangeable unit	
Electrode, Cable, and Electrolyte	
Electrode	Pt 6980
Cable	Cable: L 1 A
Burettes	T500 with 50 ml interchangeable unit T300 with 50 ml interchangeable unit T300 with 20 ml interchangeable unit
Sample Changer	TW alpha plus with sample tray TZ 1453
Stirrer	Rod Stirrer TZ 1844 with propeller blade TZ 1863
Pump	Membrane Pump MP 25
Lab Accessories	
TitriSoft, Thermometer W 5780 NN	
Laboratory Glassware: 600 ml Beakers, High Form, without spout	

A complete list of all required components incl. order numbers is available on request.





Figure 1: Complete set-up for an automated Kappa number determination



Figure 2: Titration head, view 1




Figure 3: Titration head, view 2

Reagents	
1	Deionized water
2	Potassium permanganate, 0.02 mol/L
3	Sulfuric acid, 2 mol/l
4	Potassium iodide, 1 mol/L
5	Sodium thiosulfate, 0.2 mol/l
All reagents should be in analytical grade or better.	

## Titration Performance


All instruments are controlled and the titration performed by the titration software TitrSoft. The principle structure is shown below. All steps are performed automatically by TitrSoft.

The performance of the blank determination is identical with the sample titration. The consumption at EQ is evaluated and stored as global variable.

The complete structure of the titration commands is shown in the annex of this application report.

A configured database with all required methods and settings for the automated Kappa number determination is available on request.

### Example Curve:



Dosing of 400 ml water to the sample



The sample is digested by stirring



10 minutes wait time

Simultaneous dosing of each 50 ml  
KMnO<sub>4</sub> and H<sub>2</sub>SO<sub>4</sub>



5 minutes reaction time

Temperature measurement



5 minutes reaction time

Dosing of 10 ml KI



Reaction is stopped

Titration is started



Titration is stopped at the first EQ

Calculation of the Kappa number

YSI, a Xylem brand  
1725 Brannum Lane  
Yellow Springs, OH 45387

+1.937.767.7241  
titration.yesi@xylem.com  
YSI.com



## Annex: Flow chart of all required TitrSoft titration commands

