OPERATIONS MANUAL

pH1200

LABORATORY pH/mV/Temp Instrument
English, Español, Português

a xylem brand

CONTENTS

Warranty	2
Getting started.....	3
Display.....	3
Keys.....	4
Buffer Set Selection.....	5
pH Calibration.....	6
Calibration with an ATC/Temp probe in the pH-AUTOLOCK mode.....	6
Calibration with manual temperature compensation in the pH-AUTOLOCK mode.....	7
Calibration with an ATC/Temp probe in the pH-NON-AUTOLOCK mode	8
Calibration with manual temperature compensation in the pH-NON-AUTOLOCK mode.....	9
pH Measurements	9
Measurement with an ATC/Temp probe in the pH-AUTOLOCK mode.....	9
Measurement with manual temperature compensation in the pH-AUTOLOCK mode	10
Measurement with an ATC/Temp probe in the pH-NON-AUTOLOCK mode.....	10
Measurement with manual temperature compensation in the pH-NON-AUTOLOCK mode	10
Temperature Measurements	11
mV Measurements.....	11
Measurement in the mV-AUTOLOCK mode.....	11
Measurement in the mV NON-AUTOLOCK mode	11
Troubleshooting	12
Specifications.....	13
Accessories/Spare Parts.....	14
Recycling	14
Battery Disposal.....	14
Contact Information	15
Ordering and Technical Support	15
Service Information	15

Item #601204REF
Rev A
October 2012

©2012 YSI Incc/Xylem Inc.

The YSI logo, YSI and Xylem are registered trademarks.

WARRANTY

The YSI pH1200 instrument is warranted for two (2) year from date of purchase by the end user against defects in materials and workmanship, exclusive of batteries and any damage caused by defective batteries. pH1200 cable/sensor assemblies are warranted for one (1) year from date of purchase by the end user against defects in material and workmanship. pH1200 instruments & cables are warranted for 90 days from date of purchase by the end user against defects in material and workmanship when purchased by rental agencies for rental purposes. Within the warranty period, YSI will repair or replace, at its sole discretion, free of charge, any product that YSI determines to be covered by this warranty.

To exercise this warranty, call your local YSI representative, or contact YSI Customer Service in Yellow Springs, Ohio at +1 937 767-7241, 800-897-4151 or visit www.YSI.com for a Product Return Form. Send the product and proof of purchase, transportation prepaid, to the Authorized Service Center selected by YSI. Repair or replacement will be made and the product returned, transportation prepaid. Repaired or replaced products are warranted for the balance of the original warranty period, or at least 90 days from date of repair or replacement.

LIMITATION OF WARRANTY


This Warranty does not apply to any YSI product damage or failure caused by:

- Failure to install, operate or use the product in accordance with YSI's written instructions;
- Abuse or misuse of the product;
- Failure to maintain the product in accordance with YSI's written instructions or standard industry procedure;
- Any improper repairs to the product;
- Use by you of defective or improper components or parts in servicing or repairing the product;
- Modification of the product in any way not expressly authorized by YSI.

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. YSI'S LIABILITY UNDER THIS WARRANTY IS LIMITED TO REPAIR OR REPLACEMENT OF THE PRODUCT, AND THIS SHALL BE YOUR SOLE AND EXCLUSIVE REMEDY FOR ANY DEFECTIVE PRODUCT COVERED BY THIS WARRANTY. IN NO EVENT SHALL YSI BE LIABLE FOR ANY SPECIAL, INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM ANY DEFECTIVE PRODUCT COVERED BY THIS WARRANTY.

GETTING STARTED

DISPLAY

<p>1. HOLD This will indicate that the reading is frozen during AUTOLOCK mode.</p>	<p>8. C°/F° Temperature units.</p>
<p>2. WAIT This will be displayed when the unit is still waiting for a stable reading.</p>	<p>9. AUTO AUTOLOCK mode indicator. (If the difference between a reading at the first second and tenth second is less than 0.01 pH or 0.2 mV then the instrument will "lock" the tenth reading).</p>
<p>3. LOBAT Low battery indicator. pH1200 uses 4 C-cell alkaline batteries. The battery compartment is on the bottom of the instrument and can be accessed by sliding the cover according to the direction of the arrows.</p>	<p>10. mV Measurement unit and mode indicators.</p>
<p>4. MAN Displayed if a temperature probe is not connected.</p>	<p>11. pH Measurement unit and mode indicators.</p>

<p>5. ATC ATC (Automatic Temperature Compensation) indicator will be displayed if a temperature probe is connected.</p>	<p>12. EFF. / % This will be displayed if the user is viewing the efficiency of the electrode. It is recommended to use a new electrode when the efficiency value is less than 75%.</p>
<p>6. CAL This will be displayed when the instrument enters the calibration mode.</p>	<p>13. Buffer selection This indicator will flash if the unit is not yet calibrated. This indicator will remain lit-up if the unit has been calibrated.</p>
<p>7. SECONDARY DISPLAY For temperature in C° or F° display.</p>	

KEYS

Key	Description
	<p>On/Off Press and hold this key for 5 seconds to power on or turn off the instrument. Once the unit is powered up, press the same key to turn on or off the backlight.</p>
	<p>Mode Selects display mode. Pressing this key changes the display sequentially to display pH-AUTO, mV-AUTO, pH and mV. The calibration values will not be affected by changing the display modes. During pH calibration, press Mode key to exit calibration mode.</p>

	<p>Clear (not active in mV mode) Used to clear the unit when an error signal appears. <i>It clears all calibration values stored in the internal memory.</i> Under normal use, the key will not be activated unless pressed and held for 2 seconds to prevent accidentally erasing stored memory.</p> <p>When the “Clear” key is pressed, all segments of the LCD will be on. After about 2 seconds the unit will enter the pH-AUTO mode. The “AUTO” and “CAL” will be on and one of the buffers in the pre-selected buffer set will start to flash. The unit must be calibrated again before use.</p>
	<p>Up/Down The two keys are used to manually enter the temperature values. They have no effect on the unit when operating in ATC mode.</p>
	<p>Stand/Slope The “Stand” and “Slope” keys are used for pH calibration of the unit. Pressing and holding the Stand key while turning on the power, will change the buffer set.</p>
	<p>Meas. / Effic. This key is used to bring the unit out of the AUTO condition when operating in the pH-AUTOLOCK or mV-AUTOLOCK mode. Press and hold this key for 3 seconds, the LCD will display the efficiency of the electrode.</p>

BUFFER SET SELECTION

The pH1200 has two buffer sets: 7.00, 4.01, 10.01 and 6.86, 4.00, 9.18. The meter is factory pre-set to 7.00, 4.01 and 10.01. To change the buffer set, turn off the unit, then press and hold the “Stand” key while turning on the unit again.

Note: *There is no need to repeat this procedure every time the unit is powered up unless the buffer settings need changed.*

pH CALIBRATION

The pH1200 can perform a one, two or three point calibration.

Note: *If the unit uses two or three point calibration, the first point must be 6.86/7.00. The second point can either be 4.00/4.01 or 9.18/ 10.01.*

CALIBRATION WITH AN ATC/TEMP PROBE IN THE PH-AUTOLOCK MODE

1. Turn the unit on. Press “**Clear**” key for 2 seconds, all LCD elements will light up. The meter clears all calibration values stored in internal memory.
2. Connect the pH electrode to the BNC connector and the ATC/Temp probe to the ATC/Temp connector of the unit: “ATC” will light up. “pH” and “AUTO” will light up. One of the buffers in the buffer set will flash.
3. Rinse the pH and ATC/Temp probes in distilled water and immerse them in the first buffer solution (6.86/7.00 if doing a multipoint calibration). Allow temperature reading to stabilize, then press and hold “**Stand**” key for 2 seconds to calibrate. The “WAIT” icon will flash until the unit detects a stable reading. Once the unit calibrates the first point, the selected buffer remains lit up while the remaining two buffers flash. The unit is ready to be sloped at the second buffer.

Note: *If the first buffer solution is 6.86/7.00, at this moment, press the “Mode” key. The instrument will save the calibration and exit the calibration mode. Single point calibration is complete. If the first buffer solution is 4.00/4.01, 9.18/10.00, at this moment, the instrument will save the calibration and automatically exit the calibration mode. Single point calibration is complete.*

4. Rinse the pH and ATC/Temp probe in distilled water and immerse them in the second buffer solution (either 4.00/4.01 or 9.18/10.01). Allow temperature reading to stabilize, then press “Slope” key to calibrate. The “WAIT” icon will flash until the unit detects a stable reading. Once the unit calibrates the second point, the selected two buffers light up and the remaining buffer starts to flash. The unit is ready to be sloped at the third buffer.

Note: *At this moment, press the “Mode” key and the instrument will save the calibration and exit the calibration mode. 2-point calibration is complete.*

5. Rinse the pH and ATC/Temp probe in distilled water and immerse them in the third buffer solution (either 9.18/10.01 or 4.00/4.01). Allow temperature reading to stabilize, then press "**Slope**" key to calibrate. The "WAIT" icon will flash until the unit detects a stable reading. Once the instrument calibrates the third point, the unit will save the calibration and automatically exit the calibration mode. 3-point calibration is complete.
6. The unit calculates and compensates for the pH electrode slope deviation corresponding to the values of the three calibration buffers. After calibration, press and hold "**Meas. / Effic.**" key for about 3 seconds to display the new electrode efficiency.

CALIBRATION WITH MANUAL TEMPERATURE COMPENSATION IN THE PH-AUTOLOCK MODE

1. Turn the unit on. Press "**Clear**" key for 2 seconds, all LCD elements will light up. The meter clears all calibration values stored in internal memory.
2. Connect the pH electrode to the BNC connector of the unit, "MAN" will light up. "pH" and "AUTO" will light up. One of the buffers in the buffer set will flash.
3. Rinse the pH probes in distilled water and immerse it in the first buffer solution (6.86/7.00 if doing a multipoint calibration). Adjust the temperature reading to that of the first buffer using the "**up**" or "**down**" keys (0.0 to 60.0 °C). Then press and hold "**Stand**" key for 2 seconds to calibrate. The "WAIT" icon will flash until the unit detects a stable reading. Once the unit calibrates the first point, the selected buffer remains lit while the remaining two buffers flash. The unit is ready to be sloped at the second buffer.

Note: *If the first buffer solution is 6.86/7.00, at this moment, Press the "Mode" key, the unit will save the calibration and exit calibration mode. Single point calibration is complete. If the first buffer solution is 4.00/4.01, 9.18/10.00, at this moment, the unit will automatically exit the calibration mode. Single point calibration is complete.*

4. Repeat steps 4-6 of "**Calibration with an ATC/Temp probe in the pH-AUTOLOCK mode**" for dual and three point calibration.

CALIBRATION WITH AN ATC/TEMP PROBE IN THE PH-NON-AUTOLOCK MODE

1. Turn the unit on. Press "**Clear**" key for 2 seconds, all LCD elements will light up. The meter clears all calibration values stored in internal memory.
2. Connect the pH electrode to the BNC connector and the ATC/Temp probe to the ATC/Temp connector of the unit: "ATC" will light up. Press "**Mode**" key to select "pH" mode. "pH" is on. One of the buffers in the buffer set will flash.
3. Rinse the pH and ATC/Temp probes in distilled water and immerse them in the first buffer solution (6.86/7.00 if doing a multipoint calibration). Allow temperature and pH readings to stabilize, then press and hold "**Stand**" key for 2 seconds to calibrate, the unit immediately calibrates the first point, the selected buffer remains lit while the remaining two buffers flash. The unit is ready to be sloped at the second buffer.

Note: *If the first buffer solution is 6.86/7.00, at this moment, press the "Mode" key, the unit will save the calibration and exit the calibration. Single point calibration is complete. If the first buffer solution is 4.00/4.01, 9.18/10.00, at this moment, the unit will automatically exit the calibration mode. Single point calibration is complete.*

4. Rinse the pH and ATC/Temp probe in distilled water and immerse them in the second buffer solution (either 4.00/4.01 or 9.18/10.01). Allow temperature and pH readings to stabilize, then press "**Slope**" key to calibrate. The unit immediately calibrates the second point, the other buffer solution will flash. The unit is ready to be sloped at the third buffer.

Note: *At this moment, press the "Mode" key, the unit will save the calibration and exit the calibration. Dual point calibration is complete.*

5. Rinse the pH and ATC/Temp probe in distilled water and immerse them in the third buffer solution (either 9.18/10.01 or 4.00/4.01). Allow temperature and pH readings to stabilize, then press "**Slope**" key to calibrate. The unit immediately calibrates the third point and the unit will automatically save and exit the calibration mode. Three point calibration is complete.
6. The unit calculates and compensates for the pH electrode slope deviation corresponding to the values of the three calibration buffers. After calibration, press and hold "**Meas. / Effic.**" key for about 3 seconds to display the new electrode efficiency.

CALIBRATION WITH MANUAL TEMPERATURE COMPENSATION IN THE PH-NON-AUTOLOCK MODE

1. Turn the unit on. Press **"Clear"** key for 2 seconds, all LCD elements will lit up. The meter clears all calibration values stored in internal memory.
2. Connect the pH electrode to the BNC connector of the unit, **"MAN"** will light up. Press **"Mode"** key to select **"pH"** mode. **"pH"** is on. One of the buffers in the buffer set will flash.
3. Rinse the pH probes in distilled water and immerse it in the first buffer solution (6.86/7.00 if doing a multipoint calibration). Adjust the temperature reading to that of the first buffer using the **"up"** or **"down"** keys (0.0 to 60.0 °C). Then press and hold **"Stand"** key for 2 seconds to calibrate. The unit immediately calibrates the first point, the selected buffer remains lit while the remaining two buffers flash. The unit is ready to be sloped at the second buffer.

Note: *If the first buffer solution is 6.86/7.00, at this moment, Press the "Mode" key, the unit will save the calibration and exit the calibration mode. Single point calibration is complete. If the first buffer solution is 4.00/4.01, 9.18/10.00, at this moment, the unit will immediately save and exit the calibration mode. Single point calibration is complete.*

4. Repeat steps 4-6 of **"Calibration with an ATC/Temp probe in the pH-NON- AUTOLOCK mode"** for dual and three point calibration.

pH MEASUREMENTS

To take pH measurements, the pre-selected buffer solution set must light up, indicating the unit is single-, dual-, or three-point calibrated and ready for measurements. If buffer solution set flashes, perform a pH calibration before taking measurements.

MEASUREMENT WITH AN ATC/TEMP PROBE IN THE PH-AUTOLOCK MODE

1. Connect the pH electrode to the BNC connector and the ATC/Temp probe to the ATC/Temp connector of the unit. **"ATC"** will light up.
2. Press **"Mode"** key until **"pH"** and **"AUTO"** light up.
3. Rinse the pH electrode and ATC/temp probe with distilled water and immerse in the sample. Remove air bubbles by shaking/stirring the probe.

4. Press the **"Meas."** key. **"WAIT"** will flash. The unit is waiting for a stable reading. The display will track the pH value as sensed by the pH electrode and the ATC/Temp probe.
5. When **"WAIT"** disappears, the reading is then in **"HOLD"** mode and will not respond to further changes from the sample. The pH value shown is the pH value of the sample at the displayed sample temperature. **Note:** *For samples that are inherently unstable, the unit will not AUTOLOCK. In this case, use the pH-NON-AUTOLOCK mode for measurements.*

MEASUREMENT WITH MANUAL TEMPERATURE COMPENSATION IN THE PH-AUTOLOCK MODE

1. Connect the pH electrode to the BNC connector of the unit. **"MAN"** will light up. Set unit to display the sample temperature by pressing the up and down keys (0.0 to 100.0 °C).
2. Repeat steps 2-5 of **"Measurement with an ATC/Temp probe in the pH- AUTOLOCK mode"**.

MEASUREMENT WITH AN ATC/TEMP PROBE IN THE PH-NON-AUTOLOCK MODE

1. Connect the pH electrode to the BNC connector and the ATC/Temp probe to the ATC/Temp connector of the unit. **"ATC"** will light up.
2. Press **"Mode"** key until **"pH"** lights up.
3. Rinse the pH electrode and ATC/temp probe with distilled water and immerse in the sample. Remove air bubbles by shaking/stirring the probe.
4. Allow sufficient time for the display to stabilize. The instrument will display the pH value of the sample at the displayed sample temperature.

MEASUREMENT WITH MANUAL TEMPERATURE COMPENSATION IN THE PH-NON-AUTOLOCK MODE

1. Connect the pH electrode to the BNC connector of the unit. **"MAN"** will light up. Set unit to display the sample temperature by pressing the **"up"** and **"down"** keys (0.0 to 100.0 °C).
2. Repeat steps 2-4 of **"Measurement with an ATC/Temp probe in the pH NON- AUTOLOCK mode"**.

TEMPERATURE MEASUREMENTS

The pH1200 can measure temperature independently with the ATC/temp probe without using the pH electrode. Place the ATC/temp probe in the sample. The unit will display the measured temperature.

mV MEASUREMENTS

MEASUREMENT IN THE MV-AUTOLOCK MODE

1. Connect the optional combination ORP electrode to the BNC connector of the unit.
2. Press **"Mode"** key until "mV" and "AUTO" light up.
3. Rinse electrode with distilled water and immerse in the sample. Remove air bubbles by shaking/stirring the probe.
4. Press the **"Meas."** key. "WAIT" will flash. The unit is waiting for a stable reading. The display will track the mV value as sensed by the ORP electrode
5. When "WAIT" disappears, the reading is then in "HOLD" mode and will not respond to further changes from the sample. The mV value is the sample reading.

Note: For samples that are inherently unstable, the unit will not AUTOLOCK. In this case, use the mV-NON-AUTOLOCK mode for measurements.

MEASUREMENT IN THE MV NON-AUTOLOCK MODE

1. Connect the optional combination ORP electrode to the BNC connector of the unit.
2. Press **"Mode"** key until "mV" lights up.
3. Rinse electrode with distilled water and immerse in the sample.
4. Allow sufficient time for the display to stabilize. The instrument will display the mV value of the sample.

TROUBLESHOOTING

Error on Main Display	Display Possible cause(s)	Corrective Action(s)
"Er1"	<ol style="list-style-type: none"> 1. "Stand" was pressed before the electrode and ATC/Temp probe settled to within +/-1.5 pH of the buffer value. 2. pH electrode offset is greater or less than +/- 1.5 pH. 3. pH electrode is faulty. 	<ol style="list-style-type: none"> 1. Press "Clear" key, allow sufficient time for the electrode and ATC/Temp probe to stabilize, press "Stand" again to start the calibration procedure. 2. Replace the buffer and/or the pH electrode. Press "Clear" key to recalibrate meter. 3. Replace electrode.
"Er2"	<ol style="list-style-type: none"> 1. "Slope" was pressed before the electrode and ATC/Temp probe settled to within 30% of the buffer value. 2. Buffer 4.00, 4.01, 9.18 and 10.01 are incorrect. 3. pH electrode slope is off by more than 30% of ideal slope. 	<ol style="list-style-type: none"> 1. Allow sufficient time for the electrode and ATC/Temp probe to stabilize, press "Slope" key to continue the calibration procedure. 2. Check if the correct buffer is used. 3. Replace the buffer and/or the pH electrode. Press "Clear" key to recalibrate meter.
"Er3"	<ol style="list-style-type: none"> 1. Temperature is out of the 0.0 to 60.0°C range. 	<ol style="list-style-type: none"> 1. Bring the buffer or sample temperature within range.
"over" / "undr"	<ol style="list-style-type: none"> 1. Measured pH is out of the 0.00 to 14.00 pH range. 2. Measured mV is out of the -1999.9 to 1999.9 mV range. 3. Measured temperature is out of the 0.0 to 100.0°C range. 	<ol style="list-style-type: none"> 1. Bring sample pH into the correct measuring range. 2. Bring sample ORP into the correct measuring range. 3. Bring sample temperature into the correct measuring range.

SPECIFICATIONS

Display	Range	Accuracy	Resolution
pH	0.00 to 14.00	+/- 0.01 (instrument only)	0.01 unit
Temperature	0.0 to 100°C	+/- 0.2°C	0.1°C
mV	-1999 to 1999 mV	+/- 0.05% +/- 1 digit	0.1 mV

Operating Temperature Range	0 to 50°C (32 to 122°F)
Operating Relative Humidity Range	up to 85%
Water resistance	IP-65 water resistant, splash resistant
Size	150 mm x 203 mm x 72 mm 5.9 in. x 8 in. x 2.8 in.
Weight with batteries	504 grams (1.1 lb)
Input impedance	$>10^{13} \Omega$
Power	Four 1.5 v C-size alkaline batteries or 115/230V AC adapter
pH temp compensation	Auto/Manual 0.0 to 100.0°C
pH buffer recognition	Automatic recognition US (4.01, 7.00, 10.01) or NIST (4.00, 6.86, 9.18)
pH offset recognition	+/- 1.50 pH
pH slope recognition	30%
Connector (pH/temp)	BNC/8 pin DIN
Autolock feature	Yes
Backlit Display	Yes

ACCESSORIES/SPARE PARTS

Item	Description
601201	pH Electrode, Lab, 1M
601202	Temperature Probe, Stainless Steel
601203	pH/ATC Temp Combination Probe, Lab, 1M
605170	Electrode stand
603824	pH buffer, assorted case of US buffers. 2 pints ea of buffer 4, 7 and 10

RECYCLING

YSI is committed to reducing the environmental footprint in the course of doing business. Even though materials reduction is the ultimate goal, we know there must be a concerted effort to responsibly deal with materials after they've served a long, productive life-cycle. YSI's recycling program ensures that old equipment is processed in an environmentally friendly way, reducing the amount of materials going to landfills.

- Printed Circuit Boards are sent to facilities that process and reclaim as much material for recycling as possible.
- Plastics enter a material recycling process and are not incinerated or sent to landfills.
- Batteries are removed and sent to battery recyclers for dedicated metals.

When the time comes for you to recycle, follow the easy steps outlined at ysi.com.

BATTERY DISPOSAL

The 1200 is powered by alkaline batteries which the user must remove and dispose of when the batteries no longer power the instrument. Disposal requirements vary by country and region, and users are expected to understand and follow the battery disposal requirements for their specific locale.

CONTACT INFORMATION

ORDERING AND TECHNICAL SUPPORT

Telephone: 800 897 4151 (USA)
+1 937 767 7241 (Globally)
Monday through Friday, 8:00 AM to 5:00 ET

Fax: +1 937 767 9353 (orders)
+1 937 767 1058 (technical support)

Email: environmental@ysi.com

Mail: YSI Incorporated
1725 Brannum Lane
Yellow Springs, OH 45387, USA

Internet: ysi.com

When placing an order please have the following available:

- 1.) YSI account number (if available)
- 2.) Name and phone number
- 3.) Purchase Order or Credit Card number
- 4.) Model Number or brief description
- 5.) Billing and shipping addresses
- 6.) Quantity

SERVICE INFORMATION

YSI has authorized service centers throughout the United States and Internationally. For the nearest service center information, please visit www.ysi.com and click 'Support' or contact YSI Technical Support directly at 800-897-4151 (+1 937-767-7241).

When returning a product for service, include the Product Return form with cleaning certification. The form must be completely filled out for a YSI Service Center to accept the instrument for service. The form may be downloaded from www.ysi.com by clicking on the 'Support'.

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©2012 YSI Inc./Xylem Inc.